

Western blotting

Western blotting transfers proteins from a gel onto a solid support such as a nitrocellulose or polyvinylidene fluoride (PVDF) membrane. Once transferred, individual proteins, protein–protein or protein–ligand interactions can be detected by probing for specific immunological epitopes or biochemical activity. The choice of the membrane and transfer conditions is largely dependent on the nature of the protein of interest and intended downstream applications.

Wet transfers (*Alternative A*) provide numerous opportunities for optimizing transfer conditions, including time, temperature, voltage, and buffer composition. When carefully tuned, they can achieve almost quantitative transfer for most proteins. Semi-dry transfers (*Alternative B*) are faster and use less buffer. However, they offer fewer opportunities for optimization. If transfer is incomplete, the blotting process cannot simply be extended since by that point, the thin buffer layer is typically exhausted.

The term ‘Western blot’ was coined by W. Neal Burnette (Burnette, 1981) in playful homage to the earlier ‘Southern’ and ‘Northern’ blots for DNA and RNA.

Risk assessment

- **Methanol is a REPRODUCTIVE TOXIN**
- Work with high voltage power sources
- ▷ Wear gloves, safety glasses, lab coat
- ☐ Collect as HAZARDOUS WASTE
- ☐ DO NOT wash into sewer



Reviewed: May 22, 2026

Procedures

» Choosing a Western transfer system

- (1.) Check the literature for transfer conditions of the protein of interest. Epitopes can be sensitive to the way proteins are transferred. Common starting points and optional adjustments:

Sample		Transfer		Transfer buffer			Transfer conditions		
Protein	State	System	Membrane	Recipe	MeOH	SDS	Voltage	Temp.	Time
General	Denat.	Wet	NC, 0.45 μm	Towbin	10%	Omit	100 V	4 $^{\circ}\text{C}$	30 min
Histone	Denat.	Wet	NC, 0.20 μm	Towbin	20%	Omit	30 V	4 $^{\circ}\text{C}$	70 min
Histone	Denat.	Wet	PVDF, 0.20 μm	CAPS	Omit	0.01%	100 V	4 $^{\circ}\text{C}$	60 min
Membrane proteins	Denat.	Wet	PVDF, 0.45 μm	Towbin	20%	0.10%	35 V	4 $^{\circ}\text{C}$	16 h
Transcription factors	Denat.	Semi-dry	PVDF, 0.20 μm	Bjerrum	10%	0.01%	15 V	25 $^{\circ}\text{C}$	45 min
Acidic protein	Native	Wet	PVDF, 0.20 μm	Dunn	5%	0.01%	40 V	4 $^{\circ}\text{C}$	60 min
Basic protein	Native	Wet	PVDF, 0.20 μm	CAPS	5%	0.01%	20 V	4 $^{\circ}\text{C}$	60 min

Hint: For most proteins, the defaults will work. For especially large, small, or hydrophobic proteins that transfer poorly under default conditions, adjust membrane pore size, SDS and methanol concentrations, or buffer pH as indicated in the table.

- (2.) Start by choosing the membrane material, considering the protein’s abundance and hydrophobicity, the detection method, and whether the blot will be reprobed.

Membrane	Cost	Binding capacity	Hydrophobicity	Luminescence	Fluorescence	Stripping
PVDF	50 $\text{€}/\text{cm}^2$	150–200 $\mu\text{g}/\text{cm}^2$	High	✓	Limited	✓
NC, 0.45 μm	35 $\text{€}/\text{cm}^2$	80–100 $\mu\text{g}/\text{cm}^2$	Low	✓	✓	Limited

Hint: PVDF binds more protein, making it useful for low-abundance targets. However, as proteins penetrate deeper into PVDF membranes, diffusion may lead to fuzzy bands, especially if the transfer isn’t tightly controlled. Soaking the gel in 30% PEG1000 (or higher) for 2 h can increase sensitivity 10- to 100-fold (Zeng et al., 1990).

Hint: While nitrocellulose is too brittle for repeated probing after stripping off the antibodies, it gives lower background in fluorescence applications, and often sharper bands.

- (3.) Once the membrane and transfer system are chosen, fine-tune the methanol and SDS concentrations based on the protein's size and state.

Transfer Membrane	System	Protein size	Denatured proteins		Native proteins	
			MeOH	SDS	MeOH	SDS
NC	Semi-dry	< 20 kDa	20%	Omit	5–10%	0.01%
		20–80 kDa	10%	0.01%	Discouraged	Discouraged
		> 80 kDa	10%	0.05%	Discouraged	Discouraged
	Wet	< 20 kDa	20%	Omit	5–10%	0.01%
		20–80 kDa	10%	0.05%	Discouraged	Discouraged
		> 80 kDa	10%	0.10%	Discouraged	Discouraged
PVDF	Semi-dry	< 20 kDa	20%	0.01%	5–10%	0.01%
		20–80 kDa	10%	0.01%	5%	0.01%
		> 80 kDa	5%	0.05%	Omit	0.05%
	Wet	< 20 kDa	20%	0.01%	5–10%	0.01%
		20–80 kDa	10%	0.05%	5%	0.05%
		> 80 kDa	5%	0.10%	Omit	0.10%

Hint: Methanol is thought to improve protein adsorption to nitrocellulose membranes and reduces swelling of polyacrylamide gels in low-ionic strength buffers. Too much methanol, however, may cause protein precipitation and impede transfer. Where possible, omit methanol to reduce hazardous waste disposal.

Hint: SDS is optional, but improves the transfer efficiency of proteins larger than 100 kDa or which have a tendency to precipitate. Small proteins or peptides, on the other hand, may completely penetrate through nitrocellulose membranes ("blow-through"). SDS should be omitted in this case.

- (4.) Finally, choose a transfer buffer with a pH above the isoelectric point (pI) of the most basic protein of interest. Towbin buffer is often a good starting point.

Transfer	Recipe	pH	System		Application
			Wet	Semi-dry	
Continuous	Towbin	8.3	✓	Not suitable	Works for most proteins (wet transfer only)
	Bjerrum	9.2	✓	✓	Versatile; preferred for native proteins
	Dunn	9.9	✓	Not suitable	Basic proteins
	CAPS	11.0	✓	✓	Basic proteins or prior to Edman sequencing
Discontinuous	Cathode buffer	9.4	Not applicable	✓	
	Anode buffer	10.4	Not applicable	✓	

This is why: Buffer pH determines protein mobility during transfer. Buffers like Towbin and Dunn are not recommended for semi-dry transfer due to their low conductivity and limited buffering capacity in thin stacks. Discontinuous transfer uses separate anode and cathode buffer to generate an ion gradient, but these are rarely customized by users.

🔗 [KS09]

» Preparing the transfer system

- ☐ Nitrocellulose membrane, pore size 0.20 μm or 0.45 μm ☐ Methanol, analytical grade
 – or: PVDF membrane, pore size 0.20 μm or 0.45 μm

- (1.) Freshly prepare the 1 × transfer buffer (or cathode and anode buffer for discontinuous transfers), supplementing methanol or SDS as needed. Bring buffers down to 4 °C.

Critical: Do not reuse transfer buffers! Gradient gels always require methanol in transfer buffers to prevent swelling into a trapezium shape during transfer. ←

- (2.) Pre-equilibrate the gel in transfer buffer (or cathode buffer) for 10–20 min before assembly.

This is why: Removing excess SDS and salts minimizes band distortion and blow-through. The gel will settle to its final size.

- (3.) Cut a piece of membrane to the size of the gel. Avoid touching with bare hands.

Hint: Use a paper trimmer to pre-cut membrane sheets that fit Mini (7.2 cm × 8.6 cm) or Midi Gels (13.3 cm × 8.7 cm).

- (4.) **Critical:** Briefly activate PVDF membranes by soaking in methanol, even for methanol-free transfers. 
- (5.) Hydrate the nitrocellulose or PVDF membrane in transfer buffer (or anode buffer) for 5 min. The membrane will sink on its own once fully soaked.

A > Wet transfer/Tank transfer

- | | |
|--|---|
| ☐ Nitrocellulose membrane, pore size 0.20 μm or 0.45 μm | ☐ Transfer tank |
| – or: PVDF membrane, pore size 0.20 μm or 0.45 μm | ☐ Ice pack |
| ☐ Blotting paper, 2 mm thick | ☐ Roller |
| ☐ Blotting cassette | ☐ Power supply |
| ☐ Foam pads | ☐ 1 $\times$ Transfer buffer |
| ☐ Stir bar | |

- (1.) Soak two sheets of blotting paper in transfer buffer until fully saturated.
- (2.) Assemble the gel sandwich in the blotting cassette, submerged in transfer buffer: 

Top

Foam pad

Blotting paper

Gel

Membrane

Blotting paper

Foam pad

Bottom

This is why: Placing the gel on top of the membrane allows you to visually check the alignment and remove air bubbles. The blotting paper protects the gel and membrane from concentrated ions near the electrodes. Too little paper results in lower transfer efficiency and leads to damaged electrode surfaces in the long run.

Hint: When transferring multiple gels in one sandwich, interleave them with pieces of water-soaked dialysis membrane.

- (3.) Use a roller to gently remove air bubbles trapped between the layers of the blot assembly. 
- (4.) Close the cassette and insert it into the transfer tank with the membrane facing the anode (often red).
- (5.) Add a stir bar and ice pack. Fill the tank with transfer buffer, place on a stir plate, and start stirring to ensure even temperature and ion distribution during transfer. 
- (6.) Connect the power leads. Run the transfer at one of the following regimes:  30–180 min 

- ☐ Run at constant voltage: 100 V for 30 min at 4 °C.

This is why: Constant voltage ensures that field strength remains constant, providing the most efficient transfer possible for tank blotting, but generates more heat. Use cooling systems or run at lower voltages to minimize band distortion.

- ☐ Run at constant current: 3 mA/cm² (Mini Gel, 150 mA; Midi Gel, 300 mA) for 180 min at 4 °C.

This is why: Running at constant current minimizes heat generation, but slows the transfer over time.

Hint: For proteins larger than 150 kDa, increase transfer time to 60 min; for proteins smaller than 30 kDa, decrease to 20 min. Transfer efficiency and quality generally increase with reduced power and increased transfer time. Use the ranges below as a starting point when switching to other transfer buffers:

Transfer Recipe	For 1 h		For 3 h		For 16 h	
	Voltage	Current	Voltage	Current	Voltage	Current
Towbin	50–100 V	200–400 mA	25–50 V	100–200 mA	25–40 V	40–80 mA
Bjerrum	50–100 V	200–400 mA	25–50 V	100–200 mA	25–40 V	40–80 mA
Dunn	40–80 V	200–500 mA	20–40 V	100–250 mA	10 V	40–80 mA

- (7.) Turn off the power supply and disconnect the electrical leads.
- (8.) **Optional:** Stain the gel with Coomassie blue  SOP0010 to verify transfer efficiency. 

B > **Semi-dry transfer**

- | | |
|--|--|
| ☐ PVDF membrane, pore size 0.20 μm or 0.45 μm | ☐ Power supply |
| – or: Nitrocellulose membrane, pore size 0.20 μm or 0.45 μm | ☐ 5 $\times$ Anode buffer, pH 10.4  R0115 |
| ☐ Blotting paper, 2 mm thick | ☐ 5 $\times$ Cathode buffer, pH 9.4  R0116 |
| ☐ Semi-dry blotter unit | ☐ 1 $\times$ Anode/Cathode buffer |
| ☐ Roller | |

- (1.) If using a discontinuous buffer system, soak one sheet of blotting paper in anode buffer and one in cathode buffer. For continuous systems, soak both sheets in transfer buffer.

Hint: In most cases, a continuous buffer system like Bjerrum or CAPS is preferred. Discontinuous systems are used for specific commercial protocols and require exact matching of buffer formulations.

- (2.) Assemble the transfer stack on the anode plate (typically the bottom electrode):

Top	
Blotting paper	Cathode buffer
Gel	Cathode buffer
Membrane	Anode buffer
Blotting paper	Anode buffer
Bottom	

Hint: Make sure the membrane is between the gel and the anode side so proteins migrate toward it. The positive electrode (anode) attracts the negatively charged proteins.

- (3.) Use a roller to gently remove any air bubbles trapped between the layers.

Quality assurance: Bubbles are common between the gel and membrane—press gently but thoroughly. Trapped air can lead to white spots (no transfer) on the final blot.

- (4.) Close the semi-dry blotter and connect the electrical leads. Start the transfer according to the manufacturer's instructions.

- Run at constant current: 1.5 mA/cm² (Mini Gel, 75 mA) for 15–60 min at room temperature.

 15–60 min

- (5.) Turn off the power supply and disconnect the electrical leads.

Reversible staining with Ponceau S for total protein normalization

- | | |
|--|---|
| ☐ Incubation tray | ☐ Forceps |
| – or: Incubation box | ☐ 0.01% Ponceau S solution  R0117 |

- Ponceau S is the standard for reversible membrane staining
- Amido black is sometimes reversible under harsh destaining conditions

- (1.) Carefully remove the membrane from the gel sandwich and place it into a clean incubation tray or box filled with deionized water.

Critical: Only handle membranes by the edges with clean forceps. Be careful not to touch the membrane with hands or gloves.

- (2.) Decant the water. Add Ponceau S solution until the blot is submerged. Shake for 2–10 min on an orbital shaker (100 rpm) at room temperature.

- (3.) Rinse with deionized water to remove the background stain.

- (4.) Image the membrane. Quantify total transferred protein for normalization if needed.

This is why: Total protein staining with Ponceau S provides a complementary loading and transfer control compared to single housekeeping gene antibodies, which can vary with experimental conditions. Quantify lane intensities to normalize target protein signals across samples.

- (5.) Destain in aqueous buffers or deionized water. Proceed with blocking or detection.

Storing and reusing membranes

- ☐ Plastic sheet protector
– or: Resealable plastic bag

- (1.) For short-term storage between probing rounds, keep the membrane in TBST or PBS-T at 4 °C. Wrap the container in plastic film to prevent evaporation.

Quality assurance: Membranes stored wet should be used within one to two weeks. Extended wet storage increases the risk of microbial contamination and signal degradation.



- (2.) Rinse the membrane with deionized water and place it protein side up on a clean sheet of blotting paper. Do not cover the membrane during drying.

Quality assurance: Drying preserves protein–membrane interactions and reduces risk of microbial growth. Membranes can be stored dry at 4 °C for months or even years.



- (3.) Once dry, place the membrane in a resealable plastic bag or a clean document protector. Label with date and experiment details.

- (4.) To reuse, rehydrate the membrane in water for 2 min, then transfer to TBST or another buffer.

Critical: For PVDF membranes, pre-wet in 100% methanol for 10–30 s before water to ensure proper rehydration.



Materials

Preparing the transfer system

Membrane

- Nitrocellulose membrane, pore size 0.20 µm or 0.45 µm
- or PVDF membrane, pore size 0.20 µm or 0.45 µm

Methanol, CAS 67-56-1, analytical grade

Wet transfer/Tank transfer

Membrane

- Nitrocellulose membrane, pore size 0.20 µm or 0.45 µm
- or PVDF membrane, pore size 0.20 µm or 0.45 µm

Blotting paper, 2 mm thick

Blotting cassette

Foam pads

Stir bar

Transfer tank

- e.g. Mini Trans-Blot® Cell, Bio-Rad

Ice pack

Roller

Power supply

- e.g. PowerPac™ HC, Bio-Rad

Transfer buffer, 1 ×

Semi-dry transfer

Membrane

- PVDF membrane, pore size 0.20 µm or 0.45 µm
- or Nitrocellulose membrane, pore size 0.20 µm or 0.45 µm

Blotting paper, 2 mm thick

Semi-dry blotter unit

- e.g. Trans-Blot® SD Semi-Dry Transfer Cell, Bio-Rad
- e.g. iBlot™ 2 Gel Transfer Device, Invitrogen

Roller**Power supply**

- e.g. PowerPac™ HC, Bio-Rad

Anode buffer ⚡ R0115, 5 ×, pH 10.4

Cathode buffer ⚡ R0116, 5 ×, pH 9.4

Anode/Cathode buffer, 1 ×

Reversible staining with Ponceau S for total protein normalization**Incubation vessel**

- Incubation tray
- or Incubation box

Forceps

Ponceau S solution ⚡ R0117, 0.01%

Storing and reusing membranes**Membrane storage container**

- Plastic sheet protector
- or Resealable plastic bag

Recipes**Tris-Gly transfer buffer (Towbin), pH 8.3, 10 ×**

Amount	Ingredient		Stock	Final
30.3 g	Tris base	[77-86-1]	121.14 g/mol	0.25 M
144.1 g	Glycine	[56-40-6]	75.07 g/mol	1.92 M
To 1 L	Water, reagent-grade			

Do not adjust pH.

10 × Tris-Gly transfer buffer (Towbin)

25 mM Tris base, 192 mM Glycine, pH 8.3

No GHS hazards at the indicated concentrations

- ☐ Dilute below the discharge limit
- ☐ Flush to drain with excess water

Date: Sign: R0112

Tris-Gly transfer buffer (Bjerrum), pH 9.2, 10 ×

Amount	Ingredient		Stock	Final
58.2 g	Tris base	[77-86-1]	121.14 g/mol	0.48 M
29.3 g	Glycine	[56-40-6]	75.07 g/mol	0.39 M
To 1 L	Water, reagent-grade			

Do not adjust pH. **Note:** Use for native proteins.

10 × Tris-Gly transfer buffer (Bjerrum)

48 mM Tris base, 39 mM Glycine, pH 9.2

No GHS hazards at the indicated concentrations

- ☐ Dilute below the discharge limit
- ☐ Flush to drain with excess water

Date: Sign: R0113

CAPS buffer, pH 11.0, 100 ×

Amount	Ingredient	Stock	Final
221 g	3-Cyclohexylamino-1-propanesulfonic acid (CAPS) [1135-40-6]	221.32 g/mol	1.0 M
To 1 L	Water, reagent-grade		

Adjust pH with sodium hydroxide. **Note:** Use for basic proteins. **This is why:** Some protocols use 100 mM CAPS for protein crystallography or general biochemistry, but this is not typical for Westerns.

100 × CAPS buffer

pH 11.0

**WARNING**

Skin and serious eye irritation

- ☐ Neutralise to pH 6–9
- ☐ Flush to drain with excess water

Date: Sign: R0114

Transfer buffer, 1 ×

Amount	Ingredient	Stock	Final
100 mL	Transfer buffer	10 ×	1 ×
100–200 mL	Methanol, analytical grade [67-56-1]	☐ 32.04 g/mol	10–20%
	Flammable liquid (H225); Organ toxicity, single exposure (H370, H371)		
2.5 mL	SDS ☐ R0047	20%	0.05%
	Skin irritation (H315); Serious eye damage (H318)		
To 1 L	Water, reagent-grade		

Prepare freshly. Methanol and SDS are optional. **Critical:** Methanol must be analytical grade since metallic contaminants will plate on the electrodes and increase the conductivity of the transfer buffer.

1 × Transfer buffer

1 × Transfer buffer, ☐ 10–20% Methanol,
☐ 0.05% SDS

Contains a component outside the substance collection; classify it from its SDS before use.

- ☐ Collect as ORGANIC SOLVENT WASTE
- ☐ Keep away from oxidizers

Date: Sign:

Anode buffer, pH 10.4, 5 ×

Amount	Ingredient	Stock	Final
15.1 g	Tris base [77-86-1]	121.14 g/mol	125 mM
To 1 L	Water, reagent-grade		

Adjust pH with sodium hydroxide.

5 × Anode buffer

pH 10.4

**WARNING**

Skin and serious eye irritation

- ☐ Neutralise to pH 6–9
- ☐ Dilute below the discharge limit
- ☐ Flush to drain with excess water

Date: Sign: R0115

Cathode buffer, pH 9.4, 5 ×

Amount	Ingredient	Stock	Final
15.1 g	Tris base [77-86-1]	121.14 g/mol	125 mM
26.2 g	6-Amino- <i>N</i> -caproic acid [60-32-2]	131.17 g/mol	
To 1 L	Water, reagent-grade		

Adjust pH with sodium hydroxide.

5 × Cathode buffer

pH 9.4

No GHS hazards at the indicated concentrations

- ☐ Collect as HAZARDOUS WASTE
- ☐ DO NOT wash into sewer

Date: Sign: R0116

Anode/Cathode buffer, 1 ×

Amount	Ingredient	Stock	Final
20 mL	Anode/Cathode buffer	5 ×	1 ×
10–20 mL	Methanol, analytical grade [67-56-1]	□ 32.04 g/mol	10–20%
	Flammable liquid (H225); Organ toxicity, single exposure (H370, H371)		
0.25 mL	SDS	□ R0047	20% 0.05%
	Skin irritation (H315); Serious eye damage (H318)		
To 100 mL	Water, reagent-grade		

Prepare freshly. Methanol and SDS are optional. **Critical:** Methanol must be analytical grade since metallic contaminants will plate on the electrodes and increase the conductivity of the transfer buffer.

1 × Anode/Cathode buffer

1 × Anode/Cathode buffer, □ 10–20% Methanol,
□ 0.05% SDS

Contains a component outside the substance collection; classify it from its SDS before use.

- Collect as ORGANIC SOLVENT WASTE
- Keep away from oxidizers

Date: Sign:

Ponceau S solution, 0.01%

Amount	Ingredient	Stock	Final
1 mg	Ponceau S	[6226-79-5] 760.6 g/mol	0.01%
1 mL	Acetic acid, glacial	[64-19-7] 1.74 M	1.0%
To 100 mL	Water, reagent-grade		

Note: The 0.01% solution is as sensitive as a 0.1% solution.

0.01% Ponceau S solution

0.01% Ponceau S, 1.0% Acetic acid

No GHS hazards at the indicated concentrations

- Hold for EHS review before disposal

Date: Sign: R0117

Troubleshooting**Wet transfer/Tank transfer***In Step 2:*

- Blurry or wavy bands or lateral smear
 - Ensure gel and membrane are evenly aligned and firmly pressed.
 - Use intact, uncompressed foam pads. Avoid reusing compressed blotting paper.

In Step 3:

- White spots or voids on Ponceau S stain
 - Remove trapped air bubbles by gently pressing the sandwich with a roller or conical tube.
 - Pre-wet blotting paper and pads thoroughly. Degas buffer before use if needed.

In Step 5:

- Uneven transfer efficiency across blot
 - Stir transfer buffer more vigorously during transfer.
 - Use fresh transfer buffer. Cool if appropriate.

In Step 6:

- Faint or absent bands
 - Check transfer time, voltage, and buffer composition. Verify methanol and SDS levels for the protein size.
- No visible transfer on membrane
 - Some membranes, especially PVDF, may have a designated active or coated side; ensure this side faces the gel during assembly.
 - Ensure membrane is between gel and anode.
 - Confirm black lead (cathode) is behind the gel and red lead (anode) is behind the membrane.
 - Verify electrode cables and apparatus function.

Semi-dry transfer*In Step 4:*

- Arcing or excessive heat during transfer
 - Check that all components are properly wetted and assembled with no gaps.

List of references

- W. Burnette, *Anal. Biochem.* **112**(2), 195—203 (1981).
C. Zeng, Y. Suzuki, and E. Alpert, *Anal. Biochem.* **189**(2), 197—201 (1990).
B.T. Kurien and R.H. Scofield, *Methods Mol. Biol.* **536** 9—22 (2009).

Change log

- | | | |
|------------|------------------------|--|
| 2020-03-22 | Michael Haugbro | Initial protocol. |
| 2023-01-13 | Benjamin C. Buchmuller | Adaptation as SOP. Added transfer-buffer selection guide from Carl Roth "Information Brochure on Transfer Buffers and Parameters"; troubleshooting informed by Bio-Rad "Western Blotting / Immunoblotting Introduction" (Bulletin 2895) and PhosphoSolutions "Western Blot Transfer Efficiency – The Good, the Bad, and the Ugly". |
| 2026-05-22 | Benjamin C. Buchmuller | Fix transcription error in Anode/Cathode buffer SDS final concentration (was 0.01%, now 0.05% — the prepared buffer at 0.25 mL of 20% SDS in 100 mL delivers 0.05%, well within the normal Towbin-style supplement range). |
| 2026-05-22 | Benjamin C. Buchmuller | Fix transcription error in CAPS buffer amount (was 22.1 g, now 221 g — recipe v1 delivered ~0.1 M instead of the declared 1 M, 10× under-prepared; corrected recipe yields the intended 1 M working stock that dilutes 100× to the 10 mM CAPS used for Western transfers). |

From the *Lab Protocols* collection by B. C. Buchmuller and contributors (2026). Licensed under Creative Commons Attribution Share Alike 4.0 International; see <https://creativecommons.org/licenses/by-sa/4.0/> for the terms.

For research use. Provided in good faith, without warranty or liability for any use or results. Users are responsible for compliance with local regulations and institutional policies. Where a brand name, manufacturer, or specific source is cited, it is given as an example and does not imply that other suitable equivalents are excluded.

Current when printed. Visit <https://benjbuch.github.io/check/> or scan the QR code to check for updates.



3c8d72c

