

>>

Casting native polyacrylamide gels

- | | |
|--|---|
| ☐ Gel casting apparatus | ☐ Gel comb, 15 wells |
| ☐ Gel cassette, 0.5–1.5 mm spacer | ☐ Isopropyl alcohol |

- (5.) Pour the resolving gel solution into the assembled gel cassette with the comb removed.

Gel fraction	Spacer	Mini Gel (7.2 cm × 8.6 cm)	Midi Gel (13.3 cm × 8.7 cm)
Resolving gel	1.0 mm	5.0 mL	9.0 mL
Stacking gel	1.0 mm	1.5 mL	2.5 mL

- (6.) Gently overlay the resolving gel with a few drops of isopropyl alcohol to create a smooth interface.

- (7.) Allow 60–120 min for the gel to polymerize. This ensures reproducible pore size.

⌚ 90 min

- (8.) Pour off the alcohol, rinse with water, dry the area above the resolving gel with a filter paper.



- (9.) Add ammonium persulfate and TEMED to the stacking gel solution, invert.

- (10.) Pour the stacking gel solution in the gel cassette. To avoid spilling, place the pipet tip at an inclined angle and dispense slowly. Fill to the rim before inserting the comb.

Hint: In high-percentage acrylamide gels (greater than 10% acrylamide) or when sample volumes are low, the resolving gel can provide sufficient focusing such that a stacking gel can be omitted.

- (11.) Store gels flat at 4 °C with the comb inserted for up to three months. Wrap in a wet paper towel and seal in a plastic bag.

>> Loading and running of native TBE acrylamide gels

- | | |
|---|--|
| ☐ Electrophoresis tank | ☐ 10 × Tris-borate-EDTA |
| ☐ Power supply | ☐ 2.5 M Sucrose |

- (1.) Assemble the electrophoresis cell.

- (2.) Place the gel in the electrophoresis tank, fill with TBE buffer from the same batch used for casting. Remove air bubbles trapped beneath the gel with a bent Pasteur pipette or syringe needle.



Critical: Always fill the reservoirs with the recommended amount of running buffer to prevent excessive heating.



- (3.) Remove comb by pulling straight up, slowly and gently. Flush out the wells with running buffer.

- (4.) Pre-run the gels at 4–6 V/cm for 20–30 min before loading the samples.

⌚ 30 min

This is why: This helps to stabilize the buffer system in low-ionic strength buffers and removes residual unpolymerized acrylamide and excess APS and TEMED which can interfere with migration. Thicker gels (1.5 mm) may benefit from a longer prerun to fully equilibrate buffer ions.

Hint: Pre-running is not necessary for Tris-Gly buffers or high-ionic strength buffers such as 1 × TBE.

- (5.) Prepare samples in 1 × sample buffer containing 10–15% sucrose and load into the wells.

Critical: Do not attempt to expel all of the sample from the loading device and complete loading quickly.



Hint: The detection limit per protein band is about 40 ng for a Coomassie Brilliant Blue R-250 stain or 0.5–1.0 ng for fluorescent dyes. For DNA, at least 1.0 ng is needed per band when stained with ethidium bromide.

- (6.) Connect the electrodes to the power device. Turn on the power and begin the electrophoresis run.

⌚ 60 min

- For the first 10 min allow the samples to stack using a reduced field strength of 5–10 V/cm gel length, typically 40–50 V.
- Continue at the recommended voltage until separation is complete.

Critical: Non-denaturing polyacrylamide gels are usually run at voltages between 1–8 V/cm gel length. If electrophoresis is carried out at higher voltage, differential heating in the center of the gel may cause bowing of bands or even melting of the protein and DNA complexes. In this case, run gels at 4 °C.



Note: Since TAE does not provide as much buffering capacity as TBE, TAE gels must be run more slowly.

Troubleshooting

Casting native polyacrylamide gels

In Step 7:

- Swirls in gel; Disturbed protein separation due to unequal or incomplete polymerization of the gel.
 - Ensure alcohol is poured off within one hour to avoid dehydration.
 - Let the resolving gel polymerize overnight at room temperature; pour off alcohol latest after one hour to avoid dehydration.
 - Let the gel polymerize overnight at room temperature for more consistent results.
 - If polymerization is too fast (less than 10 min), reduce APS and TEMED by 25%.
 - If polymerization is too slow (longer than 120 min), increase APS and TEMED concentration by 50%.
- Gel is brittle or too soft.
 - Cross-linker (bis-acrylamide) content is too high or too low. Get a new batch of acrylamide-bis acrylamide solution.

In Step 10:

- Poor well formation.
 - Degas the monomer solution to remove dissolved oxygen, which inhibits polymerization.
 - Use fresh catalyst stocks and adjust concentration to 0.06% TEMED and 0.08% APS if necessary.
- Webbing; excess acrylamide behind the comb teeth.
 - Use fresh APS and TEMED or slightly increase their concentration to improve polymerization at the interface.

Loading and running of native TBE acrylamide gels

In Step 2:

- Bands are uneven or wonky.
 - Ensure the gel is covered entirely with buffer while it is setting.
 - Pre-chill the running buffer to 4 °C or run in a cold room to prevent overheating in long runs.

Change log

2022-07-23	Benjamin C. Buchmuller	Initial version.
2025-03-18	Benjamin C. Buchmuller	Revised and expanded.
2025-03-30	Jeremy A. Owen	Optimized APS and TEMED concentration for stacking layer.

From the *Lab Protocols* collection by B. C. Buchmuller and contributors (2026). Licensed under Creative Commons Attribution Share Alike 4.0 International; see <https://creativecommons.org/licenses/by-sa/4.0/> for the terms.

For research use. Provided in good faith, without warranty or liability for any use or results. Users are responsible for compliance with local regulations and institutional policies.

Current when printed. Visit <https://benjbuch.github.io/check/> or scan the QR code to check for updates.



e3c6965