

Lentiviral generation of stable animal cell lines

Risk assessment

- **Work with lentiviral vectors (BSL-2+)**
- **Risk of exposure to replication-deficient virus through aerosols, sharps, skin or mucous membrane contact**
- **Risk of generation and exposure to replication-competent virus and insertional mutagenesis**
- **Work with human-derived material or transgenic cell lines (BSL-2)**

- ▷ New users must be trained and supervised by a designated experienced lab member for at least two procedures before working independently
- ▷ All work must be contained in a certified class II biosafety cabinet
- ▷ DO NOT use sharps or needles without written permission
- ▷ Wipe down all surfaces with freshly prepared 1% sodium hypochlorite or equivalent virucidal disinfectant (Opti-Cide® 3) before and after work; 70% ethanol is NOT effective to cleanup spills
- ▷ Transfer containers in secondary containment such as a sealed box or disinfected tray when moving between lab spaces
- ▷ Wear double gloves, splash goggles, disposable lab coat and sleeves

- ☐ DO NOT aspirate infectious waste using a vacuum line
- ☐ Collect liquid waste into a vessel with 20–30% bleach; ensure at least 60 min contact time before sink disposal; empty daily
- ☐ Discard pipette tips, plates, and any disposable items into 10% bleach or equivalent virucidal disinfectant (Vesphene® III); soak for 30 min before disposal; empty when filled to 80% capacity or after four days
- ☐ Dispose waste after inactivation as REGULATED MEDICAL WASTE



In the event of SPILL:

- ▷ Evacuate the area and allow aerosols to settle for 30 min before re-entry
- ▷ Soak thoroughly with freshly prepared 1% sodium hypochlorite or equivalent virucidal disinfectant (Opti-Cide® 3).

In the event of EXPOSURE:

- ▷ Wash affected area for 15 min
- ▷ Seek medical attention within 1 hour

Report to (all of):

- Principal Investigator/Supervisor
- Biosafety Officer

Reviewed: Jul 25, 2025

This is a bench card. Full protocol available online.

>> Production of virus particles

☐ 6-well plate	☐ DMEM + 10% FBS, without antibiotics	☐ Packaging plasmid
☐ Flip-top filter, 0.45 µm, PES	☐ Culture medium for the transduced cells	☐ Envelope plasmid
☐ Lenti-X™ 293T	☐ Transfer plasmid (with transgene)	

- (1.) On the day of transfection, replace medium with antibiotic-free DMEM + 10% FBS.

Critical: From this step forward, perform all steps in a certified class II biosafety cabinet. This includes transfection, medium changes, harvesting, and filtering virus-containing supernatant. Do not perform any viral manipulations on the open bench. Carry culture plates carefully using both hands. If moving between rooms or floors, place the cultures in a leak-proof secondary container lined with absorbent paper, and disinfect the exterior before transport.

- (2.) Transfect Lenti-X™ 293T cells at 70% confluency with the following plasmids using chemical transfection with Lipofectamine® 2000  SOP 0012. By default, use one well of a 6-well plate per transgene:

		Format	Plates (per well)		
Plasmid			6-well	12-well	24-well
Packaging plasmid	psPAX2		0.22 pmol	1.46 µg	530 ng
Envelope plasmid	pMD2.G (VSV-G)		0.12 pmol	0.43 µg	155 ng
Transfer plasmid	pLJM1, pLKO.1, etc.		0.28 pmol	1.20 µg	430 ng

Critical: This SOP describes small-scale virus production (less than 100 mL total supernatant with titers less than 1×10^9 cfu/mL). Larger volumes or higher titers may require additional IBC review or enhanced containment measures.

Lentiviral generation of stable animal cell lines

- (3.) *Optional:* After 18 h, replace the medium with 2 mL of the culture medium of the cells to transduce. ⌚ 18 h 📖
- (4.) *Optional:* Supplement the culture medium with 10 mM sodium butyrate to restrict cell growth.
- (5.) Harvest the virus-containing supernatant at 48 h, and optionally at 72 h and 96 h post transfection. Combine batches or store separately as needed. 🌙 📖
- Critical:** Handle supernatants with care; even replication-deficient vectors remain infectious until inactivated. Do not remove viral supernatant from the biosafety cabinet unless sealed and disinfected. ⬅
- (6.) *Optional:* Pre-clear the viral supernatant at $500 \times g$ for 5 min in a certified aerosol-tight centrifuge. ⬅
- Critical:** Centrifugation must be performed in sealed buckets and followed by immediate disinfection. Skip this step for small-scale preparations or if a certified aerosol-tight centrifuge is not available near the biosafety cabinet.
- (7.) Filter the virus-containing supernatant through a sterile $0.45 \mu\text{m}$ PES flip-top filter inside the biosafety cabinet to remove remaining packaging cells and debris. Do not use needles! ⬅
- Critical:** Flip-top filters reduce the risk of aerosol generation and spills, increasing user safety. If using syringe filters, do not force the plunger or push air through the filter as this will generate aerosols, or can cause the filter to pop off.
- (8.) Disinfect all tubes and external surfaces of containers with 1% sodium hypochlorite or equivalent virucidal disinfectant (Opti-Cide® 3) before removing them from the biosafety cabinet.

+ **Optional: PEG precipitation of viral supernatants**

☐ 🔗 R0195 5 × PEG virus precipitation solution, 100 mL

- (1.) Add 0.25 vol PEG virus precipitation solution per 1.0 vol of viral supernatant. ⚙
- (2.) Mix well by shaking for 1 min. Incubate with constant rocking at 4°C for at least 4 h. ⌚ 6 h 📖
- (3.) Spin at $1\,600 \times g$ for 1 h at 4°C in a certified aerosol-tight centrifuge. ⌚ 1 h ⬅
- Critical:** Centrifugation must be performed in sealed buckets and followed by immediate disinfection. Skip this step for small-scale preparations or if a certified aerosol-tight centrifuge is not available near the biosafety cabinet.
- (4.) Carefully remove the supernatant without disturbing the pellet. 📖
- (5.) Thoroughly resuspend the viral pellet in 0.05–0.10 vol PBS or serum-free culture medium without antibiotics by pipetting up and down. Incubate for 10 min. ⬅
- Critical:** Avoid generating bubbles which may inactivate the virus.
- (6.) *Optional:* Transfer to a microcentrifuge tube and spin at full speed for 3 min to remove cellular debris.

+ **Optional: Concentration of viral supernatants by centrifugation**

- (1.) If an ultracentrifuge is not available, place the viral supernatant on a 0.2 vol 50% sucrose cushion.
- (2.) Concentrate depending on the available instrumentation,
- In a high-speed centrifuge at $10\,000\text{--}20\,000 \times g$ for 240 min, or
 - In an ultracentrifuge at $90\,000 \times g$ for 90 min.
- (3.) Carefully remove the supernatant without disturbing the pellet. Resuspend as desired.

» **Storage of lentiviral particles**

- (1.) Store filtered viral supernatants at 4 °C for up to 24 h, or transfer in 1 mL aliquots to –80 °C storage.



Critical: Use screw-cap polypropylene vials for storage. DO NOT submerge the vials into liquid nitrogen! Label the secondary container clearly with “**BIOHAZARD: Infectious Material (Lentivirus)**”.



Quality assurance: Freezing and thawing reduces viral infectivity by 20%. However, storing the supernatant at 4 °C for more than two days reduces infectivity by up to 40–60%. Avoid repeated freeze–thaw cycles.



>> Transduction of host cells

☐ Culture medium for the transduced cells ☐ R0134 10 g/L Hexadimethrine bromide, 100 mL

- (1.) *Optional:* Determine the viral titer by ELISA, qPCR, or semi-quantitative lateral flow assays. 
- (2.) Seed cells to 25% confluency at the time of transduction. 
- (3.) Supplement 1 mL culture medium with 20 µL 10 g/L Polybrene® to a concentration of 0.2 g/L for a 20 × stock. 
- (4.) Prepare a dilution series of the viral supernatant in Polybrene®-supplemented medium. Apply to target cells to determine optimal transduction conditions. For a 6-well plate: 

Dilution	Neat	1:10	1:25	1:50	1:100	1:250
Viral supernatant	1 500 µL	150 µL	60 µL	30 µL	15 µL	6 µL
Culture medium	0 µL	1 350 µL	1 440 µL	1 470 µL	1 485 µL	1 494 µL
Culture medium + 0.2 g/L Polybrene®	75 µL	75 µL	75 µL	75 µL	75 µL	75 µL

- (5.) The next day, replace the medium with fresh culture medium.  

Safety: Do not aspirate viral medium using vacuum lines. Collect liquid waste into a flask containing bleach and incubate for at least 60 min before disposal!

- (6.) The next day (or the day after), split the cells 1:3 to 1:5 depending on cell type, and proceed with selection or analysis as appropriate. For antibiotic selection, refer to  SOP0019. 

Critical: Maintain cells in a subconfluent state throughout antibiotic selection to outgrow successfully transduced clones. 

Critical: For biosafety reasons, all transduced cells must be cultured in a designated incubator for at least **three to five additional passages**, typically over two to three weeks, to allow clearance of any potentially infectious particles. 

 Recipe (available online)  Troubleshooting (available online)  Notes (available online)

