

Selection and screening of stable mammalian cell lines

Risk assessment

- Work with human-derived material or transgenic cell lines (BSL-2)
- ▷ Wear gloves, safety glasses, lab coat
- Collect and dispose waste after inactivation as REGULATED MEDICAL WASTE



Reviewed: Apr 28, 2023

This is a bench card. Full protocol available online.

A > Antibiotic selection of mammalian cell lines

- Culture medium with antibiotics

- (1.) Determine the minimum effective antibiotic concentration for selection using a matrix titration format.
 - Prepare a 96-well plate with 100 μ L of culture medium in all wells.
 - Perform a row-wise dilution of the antibiotic: Add 100 μ L of the highest concentration to row A, and perform 1:2 serial dilutions down to row G. Leave row H as a no-antibiotic control.
 - In each row, seed column 1 with 100 μ L of a 4.0×10^4 mL⁻¹ cell suspension. Mix gently and transfer 100 μ L to column 2. Continue serial dilution across columns to column 12. Discard excess.
 - Incubate the plate under standard conditions.
 - Assess cell survival under the microscope after 3–7 d.
- (2.) Select cells using the lowest antibiotic concentration that resulted in complete killing of non-transgenic cells. Include a negative control (non-transfected) to verify selection stringency. Useful ranges:

Cell line	Blasticidin S	G418 (Geneticin)	Hygromycin B	Puromycin	Zeocin
Typical duration	10 days	14 days	4 days	4 days	10 days
293T	2–10 μ g/mL	Not determined	Not determined	1.0–5.0 μ g/mL	Not determined
HeLa	5–10 μ g/mL	Not determined	Not determined	1.0–2.0 μ g/mL	Not determined
MCF 10A	2–10 μ g/mL	150–400 μ g/mL	100 μ g/mL	0.5–2.5 μ g/mL	750 μ g/mL
iPS	10–20 μ g/mL	Not determined	Not determined	Not determined	Not determined

- (3.) Maintain selection by feeding cells regularly. Expand once resistant colonies are visible.
- (4.) *Optional:* Perform a second round of limiting dilution under selection to ensure the selected cell population is monoclonal.

B > Screening of mammalian cell lines

- (1.) Physically separate the population into clones using limiting dilution, pipette-based cell picking, or flow cytometry (FACS).
- (2.) Feed cells regularly and monitor for clonal outgrowth. Split as needed.

Critical: In 96-well plates, border wells lose volume faster than center wells. Top up edge wells regularly to prevent desiccation.
- (3.) Analyze each clone for expression of the transgene or phenotypic marker.

🔗 Recipe (available online) 🔗 Resources (available online) 🛠 Troubleshooting (available online) 📖 Notes (available online)

Selection and screening of stable mammalian cell lines

Lab Protocols — © B. C. Buchmuller. Licensed CC-BY-SA-4.0. Research use; provided as-is, without warranty. Current when printed. Visit <https://benjbuch.github.io/check/?q=698664f> or scan the QR code to read the full protocol or to check for updates.

