

# Whole cell lysates from tissue culture cells

## Risk assessment

- Work with human-derived material or transgenic cell lines (BSL-2)
- ▷ Wear gloves, safety glasses, lab coat
- Collect and dispose waste after inactivation as REGULATED MEDICAL WASTE



Reviewed: Mar 3, 2026

This is a bench card. Full protocol available online.

## >> Choosing a suitable lysis buffer

- (1.) Select a lysis buffer based on the target protein class and intended downstream application. The table below summarizes the solubility of representative protein classes in each buffer.  

| Protein                                                   | RIPA        | NP-40       | Other           |
|-----------------------------------------------------------|-------------|-------------|-----------------|
| Whole cell extracts                                       | Good        | Recommended |                 |
| Cytosolic proteins                                        | Good        | Good        | Tris-HCl        |
| Nuclear proteins                                          | Recommended | Limited     |                 |
| Mitochondrial proteins                                    | Recommended | Limited     |                 |
| Membrane-bound proteins                                   | Recommended | Good        |                 |
| High molecular weight proteins                            | Poor        | Poor        | Urea            |
| Structural proteins: Fibronectin, Keratin, Lamin, Tubulin | Limited     | Poor        | Tris-Triton     |
| Signaling proteins: EGFR, HSP90, c-Src                    | Poor        | Poor        |                 |
| Heterochromatin: H3K4me2                                  | Insoluble   | Insoluble   | Acid extraction |
| Transcription factors: GATA2                              | Insoluble   | Insoluble   |                 |
| Apoptotic markers: Cleaved Caspase-8                      | Insoluble   | Insoluble   |                 |

## >> Preparation of cell lysates

- Phosphate-buffered saline, pH 7.4      □  R0090 100 × Inhibitor cocktail

- (1.) Start with  $1 \times 10^6$ – $1 \times 10^9$  cells. Approximate the packed cell volume; this is 1 vol. Typically,  $1 \times 10^9$  cells will amount to 1–2 mL. 

**Critical:** Keep buffer volumes to a minimum to maintain high protein concentrations. 

- (2.) *For suspension cells:* Wash the cells twice with 5 vol ice-cold PBS containing the desired inhibitors. Collect the cells by centrifugation at  $400 \times g$  for 5 min at 4 °C after each step. Resuspend in lysis buffer.  
*For adherent cells:* Wash the cells twice with ice-cold PBS containing the desired inhibitors. Aspirate the growth medium and residual wash buffer.

Apply lysis buffer directly to the plate and scrape cells using a cold plastic cell scraper. Transfer lysate to a fresh tube.

*For tissues:* Weigh out the tissue specimen, cut into pieces on ice, transfer to a round-bottom micro-centrifuge tube. Freeze by immersing in liquid nitrogen. Homogenize in lysis buffer on ice.

- (3.) *Optional:* Sonicate the lysate twice for 10 s at 20 kHz on ice to disrupt cellular membranes, shear genomic DNA, and aid with protein solubilization. 

**Critical:** Since sonication can damage proteins, minimize pulse length and number of pulses. Apply short pulses and lower the power level to avoid frothing. 

## Whole cell lysates from tissue culture cells

### A > Protein extraction with RIPA lysis buffer

□  R0148 Radioimmunoprecipitation lysis buffer, 50 mL

- (1.) Add 50–100 vol ice-cold 1 × RIPA lysis buffer with inhibitors, about 1 mL per  $1 \times 10^7$  cells or 5 mg of tissue. Agitate the suspension for 20 min at 4 °C. ⌚ 20 min

**Critical:** Tissue samples will take 1–2 h or longer to complete lysis. Homogenize with an electric homogenizer or mortar. ←

- (2.) Remove the insoluble fraction by centrifugation at  $12\,000 \times g$  for 20 min at 4 °C. ⌚ 20 min
- (3.) Carefully collect the supernatant containing the soluble protein fraction in a new tube on ice.

### B > Protein extraction with NP-40 lysis buffer

□  R0149 NP-40 lysis buffer, 50 mL

- (1.) Add 50–100 vol ice-cold 1 × NP-40 lysis buffer with inhibitors, about 1 mL per  $1 \times 10^7$  cells or 5 mg of tissue. Agitate the suspension for 20 min at 4 °C. ⌚ 20 min

**Critical:** Tissue samples will take 1–2 h or longer to complete lysis. Homogenize with an electric homogenizer or mortar. ←

- (2.) Remove the insoluble fraction by centrifugation at  $12\,000 \times g$  for 20 min at 4 °C. ⌚ 20 min
- (3.) Carefully collect the supernatant containing the soluble protein fraction in a new tube on ice.

### C > Proteome extraction

□  R0151 Urea lysis buffer, 50 mL

- (1.) Add 5 vol ice-cold 1 × urea lysis buffer with inhibitors, about 100 µL per  $1 \times 10^7$  cells or 5 mg of tissue. Pipette the sample up and down to break up cell clumps and vortex gently.
- (2.) Boil the lysate at 95 °C for 5 min.
- (3.) *Optional:* Sonicate the lysate twice for 10 s at 20 kHz on ice to disrupt cellular membranes, shear genomic DNA, and aid with protein solubilization. 📖

**Critical:** Since sonication can damage proteins, minimize pulse length and number of pulses. Apply short pulses and lower the power level to avoid frothing. ←

- (4.) Remove the insoluble fraction by centrifugation at  $12\,000 \times g$  for 20 min at 4 °C. ⌚ 20 min
- (5.) Carefully collect the supernatant containing the soluble protein fraction in a new tube on ice.

 Recipe (available online)  Troubleshooting (available online)  Notes (available online)

Lab Protocols — © B. C. Buchmuller. Licensed CC-BY-SA-4.0. Research use; provided as-is, without warranty. Current when printed. Visit <https://benjbuch.github.io/check/?q=a298d81> or scan the QR code to read the full protocol or to check for updates.

