

In-gel protein staining with Coomassie Brilliant Blue

Risk assessment

- Work with flammable liquids
- Work with corrosive liquids
- ▷ Wear gloves, safety glasses, lab coat
- Collect used staining and destaining solutions as HAZARDOUS WASTE



Reviewed: Jan 16, 2023

This is a bench card. Full protocol available online.

>> Coomassie Brilliant Blue staining

- | | |
|--|---|
| ☐ Incubation tray | ☐ R0123 Coomassie destaining solution, 20 mL |
| ☐ Lint-free wipe | ☐ R0124 Coomassie staining solution, 20 mL |

- (1.) *Optional:* Place gel in a container with water. Microwave at 1 000 W for 1 min.
- (2.) Transfer gel into a tray with Coomassie staining solution. Rock gently at room temperature for 15 min. 15 min
- (3.) Discard staining solution. Rinse briefly with water to remove excess dye.
- (4.) Add Coomassie destaining solution. Rock at room temperature for 45 min. 45 min
- (5.) Pour off the destaining solution into a collection bottle for reuse. Replace with water and add a lint-free wipe in the tray to soak up free dye. Destain overnight until background disappears.

Recipe (available online) Troubleshooting (available online) Notes (available online)

Lab Protocols — © B. C. Buchmuller. Licensed CC-BY-SA-4.0. Research use; provided as-is, without warranty. Current when printed. Visit <https://benjbuch.github.io/check/?q=7de12ac> or scan the QR code to read the full protocol or to check for updates.

