

Generation of stable cell lines using DNA transposases

Risk assessment

- Work with human-derived material or transgenic cell lines (BSL-2)
- ▷ Wear gloves, safety glasses, lab coat
- Collect and dispose waste after inactivation as REGULATED MEDICAL WASTE



Reviewed: Nov 12, 2023

This is a bench card. Full protocol available online.

Transfection of host cells

- (1.) Seed cells to reach 70% confluence (50 000 cells per well in a 6-well plate) at the time of transfection.
- (2.) Prepare a mix of the transposase and transfer plasmids at a sub-stoichiometric ratio.



Plasmid		Per well	
Transposase plasmid	pPBase-RFP	0.05 pmol	0.32 µg
Transfer plasmid	Custom	0.20 pmol	1.00 µg

- (3.) Transfect the cells using a non-viral method such as Lipofectamine®  SOP0012.
- (4.) After 24 h, replace the medium.
- (5.) Begin selection or screening for stably transfected cells according to  SOP0019.



 Notes (available online)

Lab Protocols — © B. C. Buchmuller. Licensed CC-BY-SA-4.0. Research use; provided as-is, without warranty. Current when printed. Visit <https://benjbuch.github.io/check/?q=a3ffad2> or scan the QR code to read the full protocol or to check for updates.

