

Purification of nucleic acids from agarose gels

Risk assessment

- Transilluminators emit ultraviolet radiation, which can damage unprotected eyes and skin.
- Guanidine salts form TOXIC GASES with bleach or strong acid
- ▷ Wear gloves, safety glasses, lab coat
- DO NOT add bleach or acidic solutions directly into sample waste



Reviewed: Nov 26, 2023

This is a bench card. Full protocol available online.

>> Separation and excision of the desired DNA fragment from agarose gels

- (1.) Resolve the desired DNA fragment on a 0.6–1.2% agarose gel in TAE or lithium acetate buffer. Load at least 200 ng per lane. ⚠

Critical: To avoid carry-over of DNA, use fresh running buffer and leave empty lanes between samples. ←

Critical: Avoid TBE buffer if possible, since borate inhibits downstream enzymatic reactions. ←

- (2.) Visualize the gel using a blue light or long-wave UV source in a darkroom or gel box.

Critical: Always use long-wave UV (> 400 nm) and minimize exposure. Avoid pre-imaging gels under UV to prevent DNA damage and pyrimidine dimer formation. ←

- (3.) Excise the desired band with a clean surgical blade. Trim excess agarose to minimize the slice size. 📖
- (4.) Transfer the gel slice to a clean microcentrifuge tube.

A > Fast flush extraction

- (1.) Place the gel slice into a spin column. Centrifuge at $15\,000 \times g$ for 1 min. Collect the flow-through. 📖
- (2.) Rotate the column 180° and spin again at $15\,000 \times g$ for 1 min.

Quality assurance: Use a blue light flashlight to check for residual fluorescence. If DNA remains, add 50 μ L 5 mM Tris pH 8.0 and repeat the spin. Expect a faint greenish glow in the eluate. 💎

- (3.) Combine all flow-throughs. Discard the column.
- (4.) *Optional:* Further purify by silica spin column or ethanol precipitation 📄 SOP 0021. 📖

B > Extraction with guanidinium thiocyanate

□ R0135 Buffer QG (Extraction buffer), 0.5 mL □ Isopropyl alcohol

- (1.) Weigh the gel slice.
- (2.) Add 350 μ L Buffer QG per 100 mg gel. If the gel slice weighs less, supplement with electrophoresis buffer to 100 mg.
- (3.) Melt the agarose at 40–60 °C. Mix occasionally. Cool to room temperature.
- (4.) If the phenol red indicator turns red, add 10 μ L 3 M sodium acetate pH 5.0.

Critical: Do not skip this step—binding will fail at neutral or basic pH. ←

- (5.) Add 150 μ L isopropanol per 100 mg agarose. ⚠
- (6.) Bind DNA to silica 📄 SOP 0021. Skip Buffer PB. 📖

Purification of nucleic acids from agarose gels

🔗 [Recipe \(available online\)](#) 🛠️ [Troubleshooting \(available online\)](#) 📖 [Notes \(available online\)](#)

