

# In-gel protein staining with silver nitrate

Silver staining detects proteins after electrophoretic separation with high sensitivity down to 200 pg per band or spot. Unlike dyes, which require a minimum mass of protein to give detectable signal, silver staining amplifies even tiny amounts via autocatalytic silver reduction at protein-bound sites, particularly at thiols and other reactive side chains, which nucleate the formation of visible metallic silver.

This protocol describes the classical silver nitrate method (Chevallet et al., 2006), which uses formaldehyde-based image development in alkaline carbonate. Compared to other silver stains—such as silver-ammonia (Blum et al., 1987) or variants optimized for mass-spectrometry (Shevchenko et al., 1996)—this method is simple, robust, and inexpensive, making it particularly suited for analyzing many samples in parallel. However, the resulting signal is nonlinear and less compatible with downstream applications like mass spectrometry, immunodetection, or enzymatic activity assays. In those cases, Coomassie  SOP0010 or fluorescent stains may offer better quantification or functional recovery.

## Risk assessment

- **Formaldehyde is a KNOWN CARCINOGEN!**
  - **Silver nitrate is a STRONG OXIDANT that reacts violently with combustible and reducing materials such as alkynes, alkalis, or halides**
- ▷ Wear gloves, safety glasses, lab coat
- ☐ Collect formaldehyde solutions as HAZARDOUS WASTE
  - ☐ Collect silver containing solutions as SILVER WASTE
  - ☐ Precipitate silver nitrate solutions immediately with sodium chloride or ascorbic acid to prevent formation of explosive silver azide
  - ☐ DO NOT wash into sewer



Reviewed: Apr 5, 2025

## Procedures

### >> Protein fixation

- |                                                                     |                                              |
|---------------------------------------------------------------------|----------------------------------------------|
| <input type="checkbox"/> Fixation solution, 200 mL <i>&lt;R&gt;</i> | <input type="checkbox"/> 30% Ethanol, 400 mL |
|---------------------------------------------------------------------|----------------------------------------------|

- (1.) Fix the gel in 100–200 mL fixation solution for 30 min on a rocking table.

 30 min

*Critical:* Do not use aldehydes if mass spectrometry is planned. Fixation with methanol and acetic acid is sufficient but slightly less sensitive. 

*Hint:* Longer fixation (18 h) improves image quality. For mass spectrometry, fix for 20 min in a water bath slightly above room temperature to reduce crosslinking. Keep in mind that this increases background and may impair quantification.

*Hint:* Use at least 4- to 5-fold excess of fixation solution over gel volume. Multiple gels can be stained at once, but must float freely in the solution.

- (2.) Wash the gel twice in 30% ethanol for 10 min each.

 40 min

- (3.) Wash the gel twice in water for 10 min each.

 20 min

» **Sensitization and impregnation with silver nitrate**

- |                                                                              |                                                          |
|------------------------------------------------------------------------------|----------------------------------------------------------|
| <input type="checkbox"/> Sodium thiosulfate sensitizing solution, 200 mL (R) | <input type="checkbox"/> 0.2% Silver nitrate, 250 mL (R) |
| <input type="checkbox"/> Tetrathionate sensitizing solution, 200 mL (R)      |                                                          |

- (1.) Prepare four containers: one for the sensitization solution, two for rinsing with water, and one for the silver nitrate impregnation. Gels should move sequentially through each solution without delay.
- (2.) Sensitize by soaking one gel at a time:
  - For a fast, sensitive staining, incubate for 1 min in *sodium thiosulfate* sensitizing solution.
  - For consistent, linear staining, incubate for 45 min in *tetrathionate* sensitizing solution.
- (3.) Rinse the gel twice in water for 1 min each. Proceed immediately to silver impregnation.
- (4.) Impregnate with silver nitrate for 20–120 min. 

*Hint:* Gels may turn yellow. The duration of impregnation does not strongly affect the final image quality.

» **Image development**

- |                                                               |                                                 |
|---------------------------------------------------------------|-------------------------------------------------|
| <input type="checkbox"/> Silver nitrate developer, 200 mL (R) | <input type="checkbox"/> 2% Acetic acid, 200 mL |
|---------------------------------------------------------------|-------------------------------------------------|

- (1.) Wash one gel at a time by briefly dipping in water.

*Hint:* Prepare one box half-filled with water, one developer box per gel, and a shared stop solution box. Transfer gels along this line. Only one gel should be in the developer at a time.

- (2.) Transfer gel into developer. Shake until precipitate dissolves. 

*Critical:* A brown or gray precipitate typically appears within seconds. Shake immediately to redissolve it or risk particulate background. 

- (3.) When desired degree of staining is achieved, stop the reaction by transferring gel into 2% acetic acid. Incubate for at least 30 min.

*Quality assurance:* Intense bands or spots may appear a few minutes into development. Allow development to proceed until background is acceptable. No further changes will occur beyond 45 min. Aim to stop development when protein bands are clearly visible against a faint, uniform yellow-to-brown background. If the background darkens to a deep brown before bands are clearly resolved, the sensitization or impregnation step may need adjustment. 

*Note:* Multiple gels can be stacked in the same stop solution box at this stage.

- (4.) Wash gels at least twice for 30 min in water. 

60 min

**Storage of stained gels**

- (1.) Leave stained gels in water for up to several days.

*Critical:* For mass spectrometry, staining, washing, spot excision, and destaining should be performed on the same day. 

- (2.) For permanent storage, soak the gel in 10% glycerol for at least 20 min.

*This is why:* This prevents cracking during drying and helps preserve band integrity for documentation.

**Destaining for mass spectrometry**

- |                                                               |                                                                        |
|---------------------------------------------------------------|------------------------------------------------------------------------|
| <input type="checkbox"/> 5% Potassium ferricyanide, 50 mL (R) | <input type="checkbox"/> Silver nitrate destaining solution, 10 mL (R) |
|---------------------------------------------------------------|------------------------------------------------------------------------|

- (1.) Excise the band or spot of interest from the gel.
- (2.) Cover the gel slice with 150 µL destaining solution and incubate for 5–10 min. Remove the solution. 
- (3.) Rinse the gel slice five times with 150 µL water for 5 min each. Remove the water between rinses.

5–10 min

- (4.) Soak the gel slice in 200 mM ammonium hydrogencarbonate solution for 20 min.
- (5.) Rinse again with water as before.
- (6.) Proceed to in-gel digestion, or dry and store gel pieces at  $-20^{\circ}\text{C}$  until use.

⌚ 20 min

*Note:* If possible, proceed immediately to in-gel digestion to minimize protein oxidation.

🔗 [GWW+99]

## Troubleshooting

### Sensitization and impregnation with silver nitrate

#### In Step 4:

- Metallic silver deposits on the gel surface (silver mirror)
  - Use clean, grease-free glass plates when casting gels.
  - Avoid handling gels with bare or powdered gloves; always use powder-free gloves and minimize direct contact.

### Image development

#### In Step 2:

- Development does not take place or yields only a weak image
  - Silver solution was incorrectly diluted.
  - Potassium bicarbonate was mistakenly used instead of potassium carbonate.
  - For precious samples: stop the reaction, rinse thoroughly with water, ten times for 30 min each, and restart from the beginning.
- Strong yellow background.
  - Thiosulfate is exhausted or has decomposed. Replace with a fresh supply.
- Developer turns black; deposits fail to dissolve.
  - Likely thiosulfate exhaustion. Use fresh thiosulfate stock.

## Recipes

### Fixation solution

| Amount    | Ingredient             |           | Stock | Final |
|-----------|------------------------|-----------|-------|-------|
| 60 mL     | Ethanol, non-denatured | [64-17-5] | 100%  | 30%   |
| 20 mL     | Acetic acid, glacial   | [64-19-7] | 99.6% | 10%   |
| To 200 mL | Water, reagent-grade   |           |       |       |

#### Fixation solution

30% Ethanol, 10% Acetic acid



**WARNING**

Flammable liquid; Eye irritation; Skin irritation

☐ Collect as HAZARDOUS WASTE

Date: Sign: R0170

### Sodium thiosulfate, 10%

| Amount   | Ingredient           |             | Stock        | Final |
|----------|----------------------|-------------|--------------|-------|
| 1 g      | Sodium thiosulfate   | [7772-98-7] | 158.11 g/mol | 10%   |
| To 10 mL | Water, reagent-grade |             |              |       |

Stable for 2 weeks at  $4^{\circ}\text{C}$ . *Note:* This is a 630 mM stock.

#### 10% Sodium thiosulfate



Skin irritation

Date: Sign: R0125

### ***Sodium thiosulfate sensitizing solution***

| Amount    | Ingredient           |         | Stock | Final |
|-----------|----------------------|---------|-------|-------|
| 250 µL    | Sodium thiosulfate   | ☞ R0125 | 10%   | 0.01% |
| To 200 mL | Water, reagent-grade |         |       |       |

Sodium thiosulfate sensitizing solution

Date: Sign: R0126

### ***Tetrathionate sensitizing solution***

| Amount    | Ingredient                           |  | Stock       | Final  |
|-----------|--------------------------------------|--|-------------|--------|
| 10 g      | Potassium acetate [127-08-2]         |  | 98.15 g/mol | 500 mM |
| 605 mg    | Potassium tetrathionate [13932-13-3] |  | 302.5 g/mol | 10 mM  |
| 60 mL     | Ethanol, non-denatured [64-17-5]     |  | 100%        | 30%    |
| To 200 mL | Water, reagent-grade                 |  |             |        |

Prepare the day of use. Stable for 1 day at room temperature.

Tetrathionate sensitizing solution

500 mM Potassium acetate, 10 mM Potassium tetrathionate, 30% Ethanol



**DANGER**

**Reproductive toxicology;** Flammable liquid; Serious eye irritation; Skin irritation

☐ Collect as HAZARDOUS WASTE

Date: Sign: R0127

### ***Silver nitrate, 0.2%***

| Amount    | Ingredient                 |  | Stock        | Final |
|-----------|----------------------------|--|--------------|-------|
| 50 mg     | Silver nitrate [7761-88-8] |  | 169.87 g/mol | 0.2%  |
| To 250 mL | Water, reagent-grade       |  |              |       |

Store at room temperature, shielded from light.

0.2% Silver nitrate



**WARNING**



Hazardous to the aquatic environment

☐ Collect as SILVER WASTE

☐ DO NOT empty into drains

Date: Sign: R0128

### ***Silver nitrate developer***

| Amount | Ingredient                     |  | Stock        | Final  |
|--------|--------------------------------|--|--------------|--------|
| 30.4 g | Potassium carbonate [584-08-7] |  | 138.20 g/mol | 220 mM |
| 250 µL | Formaldehyde [50-00-0]         |  | 37%          | 0.01%  |
| 125 µL | Sodium thiosulfate ☞ R0125     |  | 10%          | 0.01%  |
| To 1 L | Water, reagent-grade           |  |              |        |

Prepare the day of use. Stable for 1 day at 4°C. **Critical:** Add formaldehyde no more than one hour before use to ensure consistent developer performance!

Silver nitrate developer

220 mM Potassium carbonate, 0.01% Formaldehyde, 0.01% Sodium thiosulfate



**WARNING**



**Carcinogenicity;** Serious eye irritation; Skin irritation

☐ Collect as HAZARDOUS WASTE

Date: Sign: R0129

### Potassium ferricyanide, 5%

| Amount   | Ingredient                          | Stock        | Final |
|----------|-------------------------------------|--------------|-------|
| 2.5 g    | Potassium ferricyanide [13746-66-2] | 329.24 g/mol | 5%    |
| To 50 mL | Water, reagent-grade                |              |       |

**Note:** This is a 150 mM stock solution.

#### 5% Potassium ferricyanide



**WARNING**

Hazardous to the aquatic environment; Serious eye irritation

☐ DO NOT empty into drains

Date: Sign: R0130

### Silver nitrate destaining solution

| Amount   | Ingredient                    | Stock | Final |
|----------|-------------------------------|-------|-------|
| 1.0 mL   | Potassium ferricyanide  R0130 | 5%    | 15 mM |
| 0.8 mL   | Sodium thiosulfate  R0125     | 10%   | 50 mM |
| To 10 mL | Water, reagent-grade          |       |       |

Prepare just before use. **Critical:** The yellowish solution is only active for 30 min; use immediately!

#### Silver nitrate destaining solution

15 mM Potassium ferricyanide, 50 mM Sodium thiosulfate

Date: Sign: R0171

### List of references

- H. Blum, H. Beier, and H.J. Gross, *Electrophoresis* **8**(2), 93-99 (1987).
- M. Chevallet, S. Luche, and T. Rabilloud, *Nat. Protoc.* **1**(4), 1852—1858 (2006).
- F. Gharahdaghi, C. Weinberg, C. Weinberg, D. Meagher, B. Imai, and S. Mische, *Electrophoresis* **20**(3), 601—605 (1999).
- A. Shevchenko, M. Wilm, O. Vorm, and M. Mann, *Anal. Chem.* **68**(5), 850—858 (1996).

### Change log

- 2006-11-22 T. Rabilloud Initial, peer-reviewed protocol article.
- 2023-01-16 Benjamin C. Buchmuller Adaptation as SOP.
- 2025-04-05 Benjamin C. Buchmuller Revised protocol sections.

**Open Protocol — Part of the *Lab Protocols* collection (2025) by B. C. Buchmuller and contributors.** This document is made available under the Creative Commons Attribution Share Alike 4.0 International License. To view a copy of this license, visit <https://creativecommons.org/licenses/by-sa/4.0/>.

For research use. Provided in good faith, without warranty or liability for any use or results. Users are responsible for compliance with local regulations and institutional policies.

Current when printed. Visit <https://benjbuch.github.io/check/> or scan the QR code to check for updates.



c29d2eb

